

Does Measurement Site Impact Lower Extremity Edema Assessment When Evaluated Using Tissue Dielectric Constant Values?

Kelvin Chan¹, Anthony Roselli², Andrea A. Astudillo³, Elham Shams⁴, Harvey Mayrovitz⁵

1. College of Allopathic Medicine, Nova Southeastern University, Davie, USA 2. College of Allopathic Medicine, Nova Southeastern University Dr. Kiran C. Patel College of Allopathic Medicine, Davie, USA 3. Osteopathic Medicine, Nova Southeastern University Dr. Kiran C. Patel College of Osteopathic Medicine, Davie, USA 4. N/A, Nova Southeastern University Dr. Kiran C. Patel College of Osteopathic Medicine, Davie, USA 5. Department of Medical Education, Division of Physiology, Nova Southeastern University Dr. Kiran C. Patel College of Allopathic C. Patel College of Allopathic Medicine, Davie, USA 5. Department of Medical Education, Division of Physiology, Nova Southeastern University Dr. Kiran C. Patel College of Allopathic Medicine, Davie, USA 5.

Corresponding author: Kelvin Chan, kc2178@mynsu.nove.edu

Categories: Cardiology, Internal Medicine, Obstetrics/Gynecology

Keywords: lymphedema measurement, edema measurement, lymphedema diagnosis, edema diagnosis, lower extremity edema, lower extremity lymphedema, lymphedema, edema

How to cite this poster

Chan K, Roselii A, Astudillo A A, et al. (2021) Does Measurement Site Impact Lower Extremity Edema Assessment When Evaluated Using Tissue Dielectric Constant Values?. Cureus 13(9): e.

Abstract

Background:

Lower extremity edema/lymphedema is the abnormal accumulation of interstitial fluid due to various etiologies. Quantification and tracking of changes are clinically useful. A novel method uses tissue dielectric constant (TDC) measurements to assess local tissue fluid by touching the skin for a few seconds. A previous study indicated that a leg/forearm TDC ratio (1.35) is a threshold for a lower extremity edema. This study's goal was to compare outcome differences when measuring medial vs. lateral leg sites for this diagnosis.

Methods:

Twelve women with a confirmed diagnosis of lower extremity edema/lymphedema were evaluated after signing an IRB approved consent. Triplicate TDC values were measured eight cm proximal to the medial malleolus of the most edematous leg on both medial and lateral sides. A site on the ipsilateral forearm was measured as a reference. Differences between medial and lateral absolute values and leg-to-arm ratios were evaluated using Wilcoxon signed-rank test with p-value <0.05 being accepted as statistically significant.

Results:

The group mean age was 76.8 \pm 12.4 years with a BMI of 31.1 \pm 10.1 kg/m². TDC values at the medial and lateral sites were not statistically different (61.8 \pm 10.4 vs. 55.8 \pm 12.6, p = 0.595). TDC leg-to-arm ratios at the medial and lateral sites were not statistically different (1.843 \pm 0.299 vs. 1.732 \pm 0.249, p = 0.585).

Discussion:

TDC is a unique way to rapidly evaluate localized edema/lymphedema. Present results indicate measurements at either medial or lateral leg sites yield clinical assessments with similar outcomes. This indicates the method can offer rapid quantitative assessments of lower extremity edema/lymphedema.

Open Access Published 09/10/2021

Copyright

© Copyright 2021 Chan et al. This is an open access poster distributed under the terms of the Creative Commons Attribution License CC-BY 4.0., which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Distributed under Creative Commons CC-BY 4.0